



Preva Radiation Dose Table

10/30/2006

Instruments Used

Radcal Model 9010	s/n	90-2215
Radcal Model 4083	s/n	43-0221
Ion Chamber Model 1C	s/n	10294

Preva Tube Head	s/n	TA60046
-----------------	-----	---------

Rectangular cone

measured kVp: measured with Radcal 4083

measured mA: read from Preva display

ion chamber placed at the end of the cone

All exposures are with 0.100 second exposure times.

kVp	mA	mR
60	4	24.6
60	5	29.6
60	6	33.8
60	7	38.1
65	4	30.7
65	5	36.2
65	6	41.8
65	7	46.3
70	4	36.2
70	5	41.9
70	6	48.3